## PRODUCT INFORMATION

### DNA Remover™

#### Decontamination reagent for the removal of DNA

##### Background

The presence of contaminant DNA at a PCR work area can result in false positive results and misinterpretation. DNA contaminations are usually sustainable and its removal extremely wasteful. DNA Remover™ is intended for use at PCR workstations. This cleansing solution is highly active against DNA contamination. The solution contains surfactant and a DNA-destroying agent. DNA Remover™ is stable and heat resistant.

DNA Remover™ is *ready-to-use* for eliminating DNA from any surface. By following the decontamination instructions below, DNA is completely destroyed and removed. DNA Remover™ should be stored at room temperature. At colder temperatures a precipitate may form which can be easily brought into solution at 37°C. Storage for 2 weeks at 65°C does not reduce the quality of the product.

##### Protection and Precaution Information

Eye contact and prolonged contact of DNA Remover™ with skin may cause irritation. Therefore, safety glasses and disposable gloves should be worn while handling the reagent.

DNA Remover™ should be applied on glass, ceramic, plastic, rubber, steel and precious metal. DNA Remover™ can not be used for the cleaning of light or non-ferrous metals. To avoid damage or discoloration, it is recommended to spot test sensitive surfaces prior to use.

##### Instructions for Use

For the decontamination of smooth, non-porous **surfaces** spray DNA Remover™ directly on the surface, soak for 1 minute and dry with a paper towel. Then rinse immediately and thoroughly with clean water and dry with a clean paper towel.

For the treatment of **laboratory devices** apply DNA Remover™ to a paper towel and wipe all exposed surfaces of the apparatus thoroughly, soak for 1 minute and wipe with a wet paper towel. To clean small parts, briefly soak them for 1 minute in DNA Remover™ rinse with clean water and dry.

Contamination of **pipettors** may even occur using filtered tips. For the decontamination follow the manufacturer's instructions and remove the shaft from the pipettor. Remove seals and gaskets from the shaft. Soak the shaft for 1 minute in DNA Remover™. Rinse immediately and thoroughly with clean water, dry and reassemble.

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Minerva Biolabs GmbH · Landhausring 7 · 12683 Berlin · Germany