

# Fluorescence *IN SITU* Hybridisation Protocol in Paraffin-Embedded Tissues Sections



## Requirements (not provided)

Reagents
Phosphate Buffered Saline concentrate tablets
Sodium Thiocyanate
Tween 20
Deionised Formamide
Pepsin
Glycine
Paraformaldehyde
Standard Saline Citrate (SSC) Stock Solution
Xylene

**Solutions to be prepared:** PBS  
Graded Alcohols  
Sodium Thiocyanate Solution  
Pepsin solution  
Glycine  
Paraformaldehyde Solution  
Formamide Wash Solution  
Stringency Wash  
Detergent Wash

**PBS:** Prepare at single and double concentration

**Graded Alcohols:** Use Analar grade 'absolute' ethanol and purified water to prepare 95%, 80% 69% and 30% alcohols.

**Sodium Thiocyanate:** Dissolve 16g sodium thiocyanate in 200ml purified water

**Pepsin solution:** Dissolve 0.8g pepsin in 200ml of 0.1M HCl just before use.

**Glycine solution:** Dissolve 0.4g glycine in in double concentration PBS

**Paraformaldehyde solution:** Dissolve 8g of paraformaldehyde (care!) in 200ml of PBS at 80°C, cool to room temp before use. Use on day of preparation.

**Formamide wash solution:** 50ml Deionised formamide mixed with 50ml 2XSSC

**Stringency wash solution:** 2XSSC. Diluted from stock SSC

**Detergent wash solution:** Add 0.1ml of 10% Tween 20 to 200ml 4XSSC.

**This product is for research use only**

Cambio Ltd, The Irwin Centre, Scotland Road, Dry Drayton, Cambridge, UK, CB23 8AR  
October 2009

Version: February 2010

# Fluorescence *IN SITU* Hybridisation Protocol in Paraffin-Embedded Tissues Sections



## Procedure: Probe preparation & denaturation

1. Collect 3-6 micron tissue sections on coated slides (such as Fisher Superfrost. Dry overnight at 37°C
2. Dewax in Xylene, 3 x 5 min each
3. Rehydrate through graded alcohols to water
4. Incubate with Sodium thiocyanate solution for 10 mins at 80°C (Care!)
5. Wash in PBS
6. Incubate in Pepsin solution for 10 min at 37°C (see notes)
7. Quench the pepsin in Glycine solution
8. Wash in PBS
9. Post-fix in paraformaldehyde solution for 2 min
10. Wash well in PBS: 3 changes over 15 min
11. Dehydrate through graded alcohols then air dry
12. Remove the pre-diluted Chromosome Paint from the freezer, mix well and warm to 37°C
13. Apply 10-15µl paint mix to the centre of the slide
14. Cover with a glass coverslip (22X40mm) and seal with Fixogum Rubber Cement
15. Denature the sealed slide at the required temperature (See below) for 10 min
16. Place the slide horizontally in a humid chamber and hybridise overnight at 37°C
17. Carefully peel of the rubber cement and remove the coverslip
18. Wash in Formamide Wash Solution at 37°C for 3 changes of 5 min each
19. Wash with Stringency Wash solution at 37°C for 3 changes over 15 min
20. Wash with Detergent Wash Solution at 37°C for 10 min
21. Wash with PBS, 3 changes over 15 min
22. Mount in reagent MD (DAPI & Mountant) and examine

**This product is for research use only**

Cambio Ltd, The Irwin Centre, Scotland Road, Dry Drayton, Cambridge, UK, CB23 8AR  
October 2009

Version: February 2010

# Fluorescence *IN SITU* Hybridisation Protocol in Paraffin-Embedded Tissues Sections



## Notes:

Many of the reagents and solutions require specific precautions: read the products inserts datasheets.

Digestion times for Pepsin solution need to be determined for individual tissues. This will depend upon the type and length of fixation as well as the tissue type itself.

We recommend that mouse paints are denatured at 60°C and Human at 80°C.

The Sodium thiocyanate step appears to be crucial. We recommend that this step is kept constant, but explore different Pepsin digestion times. Pepsin is known to autodigest. We recommend that pepsin powder is added to 0.1M HCl whilst slides are washing in PBS after step 4.

It is possible to carry out immunohistochemistry before chromosome detection, and to use indirect methods to visualise the chromosome signals [See Poulson et al. (2001). *The Journal of Pathology*. 195: 229-235].

This method was prepared Rosemary Jeffery and Richard Poulson of the Histopathology Unit, Cancer Research UK, London.

**This product is for research use only**

Cambio Ltd, The Irwin Centre, Scotland Road, Dry Drayton, Cambridge, UK, CB23 8AR  
October 2009

Version: February 2010