

# Protocol P (PRINS Probe) Fluorescence *IN SITU* Hybridisation Protocol



(1490-T)

## Approx time:

Probe Preparation: 1 hr

<b>Solutions to be prepared:</b>	Stop Buffer
----------------------------------	-------------

**Stop Buffer:** 50mM NaCl  
50mM EDTA, pH8.0

**Reaction Mix** contains dNTPs, primer, glycerol and reaction buffer.

## Procedure:

1. Slides can be made up to 24 hours in advance dehydrated through an ethanol series air dried and stored at 4°C. No pretreatment is necessary.
2. Warm Reaction mix to 37°C for 5 min and mix well.
3. Combine 25µl of reaction mix, 1nmole of Cy3 or FITC-dUTP, or 0.04mmole biotin-dUTP and 2 units of Taq polymerase per slide.

**Note:** Labels and Taq polymerase are not supplied with this product.

4. Add to slide. Seal with Fixogum and dry for 5 min in an oven at 37°C.
5. Place slide on a preheated block at 94°C for 5 min.

**Note:** Slides have to be in direct contact with the hot surface.

6. Transfer to an oven at 58°C for 30 min.
7. Wash in stop buffer for 5 min at 58°C.
8. Wash in stop buffer for 5 min at room temperature.
9. Dehydrate through an ethanol series and air dry.
10. Dilute 1µl Counterstain 1 (DAPI) with 9µl distilled water. Add 5µl diluted counterstain to 200µl Mountant. Mix well.
11. Mount slides with 20µl Mountant/Counterstain. Overlay with a coverslip and seal with nail varnish.
12. Slides should be stored in the dark at 4°C.
13. View slides using epifluorescence filters as appropriate.

**This product is for research use only**

Cambio Ltd, The Irwin Centre, Scotland Road, Dry Drayton, Cambridge, UK, CB23 8AR  
October 2009

Version: February 2010