



Simplifying Genomics



T7 R&DNA™ Polymerase

Cat. No. D7P9205K

1. Introduction

T7 R&DNA™ Polymerase* and the corresponding wild-type T7 RNA polymerase are useful for *in vitro* synthesis of defined "transcripts" that are complementary to nucleic acids cloned into a plasmid or other vector downstream from a T7 RNA polymerase promoter. In contrast to the corresponding wild-type T7 RNA polymerase, the T7 R&DNA Polymerase can incorporate 2'-deoxyribonucleoside-5'-triphosphates or other 2'-modified triphosphates such as ribonuclease-resistant 2'-fluoro-ribonucleoside-5'-triphosphates, in addition to the canonical ribonucleoside-5'-triphosphates. The ability of this mutant enzyme to incorporate various non-canonical 2'-ribonucleotides permits either primed or unprimed *in vitro* synthesis of "transcripts" composed of rNMPs, dNMPs, modified 2'-NMPs, or of mixed dNMP/rNMP, or 2'-modified-NMP/rNMP composition for a variety of applications.

2. Product Designations and Kit Components

Product	Kit Size	Catalog Number	Reagent Description	Part Numbers	Volume
T7 R&DNA™ Polymerase	5,000 Units	D7P9205K	T7 R&DNA™ Polymerase (50 U/ μ L)	E0084-50D1	100 μ L
			DTT (100 mM, dithiothreitol)	SS000065-D7	400 μ L
			R&DNA™ Polymerase 5X Reaction Buffer	SS000333-D1	1 mL

3. Product Specifications

Storage: Store only at -20°C in a freezer without a defrost cycle.

Storage Buffer: T7 R&DNA Polymerase is supplied in a 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 100 mM NaCl, 1.0 mM DTT, 0.1 mM EDTA, and 0.1% Triton® X-100.

Unit Definition: One unit converts 1 nanomole of ribonucleoside triphosphate (NTP) into acid-insoluble material in 60 minutes at 37°C.

R&DNA Polymerase 5X Reaction Buffer: 0.2 M Tris-HCl (pH 7.5), 30 mM MgCl₂, 50 mM NaCl, and 10 mM spermidine. DTT and NTPs must also be added to the final reaction.

Quality Control: T7 R&DNA Polymerase is function-tested for RNA polymerase activity and incorporation of dCTP in two independent reactions containing 40 mM Tris-HCl (pH 7.5), 6 mM MgCl₂, 10 mM NaCl, 2 mM spermidine, 10 mM DTT, 1 or 5 μ g linearized DNA template, 0.5 mM of each NTP or 0.5 mM of each ATP/GTP/UTP/dCTP and varying amounts of enzyme.

Contaminating Activity Assays: T7 R&DNA Polymerase is free of detectable exo- and endonuclease, RNase, and *E. coli* RNA polymerase activities.

4. Suggested Reaction Protocol

T7 R&DNA Polymerase will incorporate many different 2'-deoxy NTPs and 2'-modified NTPs (except for 2'-O-methyl UTP or CTP). The efficiency of incorporation, yield, and optimal reaction conditions are likely vary depending on the specific NTP substitution.

In general,

- A. 2'-fluoro-dNTPs (2'-Fl-dNTPs) are incorporated more efficiently and produce higher yields of full-length chimeric RNA/DNA than the 2'-dNTPs.
- B. Complete substitution of one 2'-Fl-dNTP (or one dNTP) for a rNTP in a T7 R&DNA Polymerase reaction results in a slight decrease in yield. Additional substitutions of NTPs will subsequently reduce yields of transcript produced. Substitution of all four NTPs will result in extremely low yields of transcript and is not recommended.

Protocol

1. Combine the following reaction components on ice in the order given:

x μ L	RNase-Free water
1 μ g	linearized template DNA with appropriate promoter
4 μ L	R&DNA Polymerase 5X Reaction Buffer
1 μ L	of each 10 mM NTP or 2'-Fl-dNTP or dNTP or other modified-NTP or dNTP
2 μ L	100 mM DTT
0.5 μ L	T7 R&DNA Polymerase
<hr/>	
20 μ L	Total reaction volume

2. Incubate at 37°C for 2 hours.
3. Isolate or purify the reaction products by method of choice for use in subsequent experiments.

5. Related References

1. Sousa, R. and Padilla, R. (1995) *EMBO J.* **14** (18), 4609.
2. Padilla, R. and Sousa, R. (1999) *Nucl. Acids Res.* **27** (6), 1561.
3. Huang, Y. et al., (1997) *Biochemistry* **36**, 8231.

*Use of T7 R&DNA Polymerase, DuraScribe® T7 Transcription Kits are covered by U.S. Patents 5,849,546; 6,107,037; 6,596,494; 7,291,487; or 7,309,570; and other issued or pending patent applications assigned to or exclusively licensed to Epicentre® (an Illumina® company). Purchase of a product from Lucigen that is covered by a claim of these patents is accompanied by a limited non-exclusive license for a purchaser to use the purchased product to synthesize nucleic acids with non-canonical bases solely for non-commercial research. A particular use is considered non-commercial only if said use is not used, directly or indirectly, for the identification, development, or production of a therapeutic, diagnostic or other commercial product. A for-profit organization requires a license from Epicentre® (an Illumina® company) in order to use a product covered by a claim of these patents to identify, develop or make a therapeutic, diagnostic or other commercial product, and a not-for-profit organization requires a license if the product, the research, or the result of screening or research is transferred to or obtained for or on behalf of a for-profit organization. Please contact Lucigen for information related to licenses.

Epicentre is a trademark of Illumina, Inc. and/or its affiliate(s) in the U.S. and other countries, and is used under license.

DuraScribe is a registered trademark of Lucigen.

Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.