

ATW0056

## APTAMER INFORMATION

### CD63 A4WWZ

#### Description:

- Identifiers: CD63 A4WWZ
- Number of DNA nucleotides: 32 bases
- Target for selection: Recombinant CD63 human protein (MyBioSource Cat# 968215)


Aptamer was selected at Base Pair Biotechnologies from a randomized 32-mer DNA library against human CD63 protein. Proprietary methods were then used to select the aptamer.

#### Aptamer folding instruction before use:

Once the aptamer is in its working concentration, it needs to be heated to 85-95 °C for 5 minutes, and then cooled to room temperature (~ 15 minutes) before use.

#### Validation data with CD63 protein by Microscale Thermophoresis (MST) method:

- Buffer used for validation: 1X PBS, pH = 7.4, 1 mM MgCl<sub>2</sub>, 0.05% Tween
- Average  $K_d$ : 40.0 ± 2.1 nM



**Figure 1.** Aptamer-CD63 binding. The fraction of bound aptamer is plotted versus the titrated target concentration,  $R^2 = 0.995$ .